



### MODEL K9901C

# HIGH FREQUENCY SHOCK TUBE

- Measures resonant frequency of pressure sensors
- Estimates high frequency response of pressure sensors
- Enables high frequency sensitivity calculation

### **TYPICAL APPLICATIONS**

- Explosive ordnance disposal suit research
- Body armor research
- Research & development




## LAB ACCURACY TO THE FIELD

The Model K9901C Shock Tube performs dynamic characterizations on a wide range of pressure sensors. The Model K9901C produces a shock wave by the sudden burst of a replaceable diaphragm. By measuring the timing of the shock wave, as well as measuring the ambient temperature, pressure, and relative humidity, engineers can calculate the sensitivity of a test sensor mounted in the end of the shock tube. The frequency response and resonant frequency of the test sensor may also be estimated as a result of this shock wave. In addition to the shock tube hardware, K9901C includes a pre-installed data acquisition system and software.

This high speed data acquisition software:

- Records, saves, and recalls pressure and weather data
- Maintains a database of equipment used for each test
- Controls data acquisition hardware parameters

APPROXIMATE SHOCK PRESSURES AS A FUNCTION OF DIAPHRAGM						
Required Driver Pressure (psi [kPa])	Incident Shock Pressure (psi [kPa])	Reflected Overpressure (psi [kPa])	Diaphragm Material	Diaphragm Material Thickness (in [mm])	Driver Gas	
1400 [9700]	180 [1200]	1050 [7200]	Aluminum	0.064 [1.6]	Helium	
880 [6100]	144 [1000]	790 [5400]	Aluminum	0.04 [1]	Helium	
760 [5200]	130 [900]	690 [4800]	Aluminum	0.032 [0.8]	Helium	
423 [2900]	90 [600]	435 [3000]	Aluminum	0.02 [0.5]	Helium	
760 [5200]	58 [400]	245 [1700]	Aluminum	0.032 [0.8]	Nitrogen	
225 [1600]	32 [200]	110 [800]	Mylar [1]	0.01 [0.3]	Nitrogen	
7 [50]	3 [20]	6.6 [50]	Mylar [1]	0.001 [0.03]	Air	



#### **Shock Tube Diagram**

SPECIFICATIONS						
Physical						
Driver Section Dimensions (Inner Diameter x Length)	2 in x 12 in	50 mm x 305 mm				
Driven Section Dimensions (Inner Diameter x Length)	2 in x 72 in	50 mm x 1 829 mm				
Driven Gas	Dry air					
Driver Gas	Dry air, helium, other inert gas					
Mounting Stand Footprint (Length x Width)	96 in x 11 in	2 438 mm x 280 mm				
Incident Position Pressure Ports	QTY. 2					
Reflected Position Pressure Ports	QTY. 1–2 <sup>[2]</sup>					

Mylar diaphragms may be stacked for higher pressure levels.
Two separate end plates are included. One with a single port, and another with dual ports for side-by-side comparison testing.

#### **SUPPLIED ACCESSORIES**

Data Acquisition System and Pre-installed Software

ICP® Pressure Sensors, Signal Conditioning, and Cabling

Mounting Adaptors for Sensors with PCB 112 Series and 109 Series Form Factor

Consumable Diaphragms



### 10310 Aerohub Boulevard, Cincinnati, OH 45215 USA

modalshop.com | info@modalshop.com | 800 860 4867 | +1 513 351 9919

© 2021 PCB Piezotronics - all rights reserved. PCB Piezotronics is a wholly-owned subsidiary of Amphenol Corporation. Endevco is an assumed name of PCB Piezotronics of North Carolina, Inc., which is a wholly-owned subsidiary of PCB Piezotronics, Inc. Accumetrics, Inc. and The Modal Shop, Inc. are wholly-owned subsidiaries of PCB Piezotronics, Inc. The Modal Shop, Inc. are wholly-owned subsidiaries of PCB Piezotronics, Inc. The Modal Shop, Inc. are wholly-owned subsidiaries of PCB Piezotronics, Inc. Intervention is provided herein, the company names and product names used in this document may be the registered trademarks or unregistered trademarks of PCB Piezotronics, Inc., PCB Piezotronics of North Carolina, Inc. (d/b/a Endevco), The Modal Shop, Inc. or Accumetrics, Inc. Detailed trademark ownership information is available at www.pcb.com/trademarkownership.